

Product information

FD 366 sensitive wipes

disinfection of sensitive surfaces



FD 366 sensitive wipes at a glance

- Practical disinfecting wipes for fast-acting disinfection of sensitive surfaces on invasive and non-invasive medical devices in hospitals, private practices and laboratories.
- Above-average material compatibility. Particularly suitable for sensitive surfaces, such as acrylic glass, plastics, and synthetic leather coverings on treatment units.
- In accordance with the VAH, the exposure time is 1 min. (bactericidal / yeasticidal).
- Broad spectrum of activity: bactericidal, tuberculocidal, yeasticidal, virucidal activity against enveloped viruses (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses.
- Tested in accordance with current VAH methods and European standards.
- Tested in 4-field test in accordance with VAH and EN 16615.
- FD 366 sensitive is included in the VAH list and IHO list of disinfectants.
- The active substance solution has been tested under low and high load in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 16615, as well as DVV/RKI guidelines.
- Wipe size: 200 x 300 mm.

Properties

FD 366 sensitive wipes from Dürr System Hygiene are aldehyde-free, alcohol-based ready-to-use wipes for quick-acting disinfection of sensitive surfaces on invasive and non-invasive medical devices and medical furnishings in private practices, laboratories and hospitals. FD 366 sensitive wipes are distinguished by fast efficiency and excellent material compatibility. Disinfected surfaces dry rapidly.

Product composition

FD 366 sensitive wipes contain FD 366 sensitive, a quick-acting combination of alcohols and quaternary ammonium compounds as well as excipients. 100 g active substance solution contain 17 g 1-propanol, 0.125 g alkyl-benzyl-dimethylammonium chloride, excipients and water.

Antimicrobial effect

FD 366 sensitive wipes have bactericidal, tuberculocidal, yeasticidal, virucidal activity against enveloped viruses (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses. FD 366 sensitive is included in the VAH list and IHO list of disinfectants. The active substance solution has been tested under low and high load in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 16615, as well as DVV/RKI guidelines. FD 366 sensitive wipes are tested as ready-to-use wipes in 4-field test in accordance with VAH or EN 16615.

Directions for use

Remove FD 366 sensitive wipes individually. Wipe the surfaces and objects to be disinfected and allow the product to take effect. Disinfect visibly contaminated surfaces with FD 366 sensitive wipes and remove contamination. Disinfect cleaned surfaces again. Make sure that wetting is complete throughout the entire action time. Close the wipe dispenser immediately following use. In accordance with VAH, the action time is 1 minute (bactericidal / yeasticidal), 2 minutes for Tb bacteria. Inactivation of viruses takes place within 30 seconds to 5 minutes.

Environmental impact

The packaging consists of a composite material and can be thermally recycled. Dispose of used wipes with the domestic waste. Do not flush.

For further information see safety data sheet.

Shelf life

Do not use the product after the expiry date printed on the packaging.

After opening, 3 months

Product information

FD 366 sensitive wipes

disinfection of sensitive surfaces

Container sizes

Dispenser box with 100 wipes
4 refill packs with 100 wipes

Storage

Store in a well-ventilated place. Keep container tightly closed.

General notes

Disinfect the exterior surfaces of medical devices, such as dental curing lights, nozzles, handpieces, and contra-angle handpieces, prior to processing as required in accordance with the manufacturers' instructions. In the case of sensitive surfaces, test for compatibility in an inconspicuous place. Changes in product colour/odour may occur, especially when stored in sunlight. However, these changes in colour and/or odour have no impact on the disinfection efficacy of the product. Use exclusively by qualified personnel. All serious events that occur in association with the product must be immediately reported to the manufacturer and the competent authority of your member state.

Manual processing of the dispenser box and lid: Mechanically scrub and remove any deposits or residues. Thoroughly rinse off

dispenser box including lid inside and out with hot tap water (> 55 °C – NB: danger of scalding!) and dry with disposable wipe. Then, using a disposable wipe with alcoholic fast-acting disinfectant, completely wet and thoroughly disinfect dispenser box and lid (allow sufficient time of action!). Make sure to dry the dispenser and lid completely.

Alternative automated processing: Process and dry the dispenser and lid in a chemothermal cleaning and disinfection system at up to 93 °C.

Hazard warnings

FD 366 sensitive wipes are classified and labelled according to the CLP Regulation. For further information see safety data sheet.

Flash point: 33 °C

Independent expert opinions – in-house investigations

The expert opinions are available upon request.



Application	Contact time
Recommendation for use for surface disinfection in accordance with VAH	Bactericidal / yeasticidal 1 min.

Type of activity	Tested according to	Contact time
Bacteria and yeasts		
Bactericidal / yeasticidal	EN 13727, EN 13624, EN 16615	1 min.
Tuberculocidal	EN 14348	2 min.
Viruses		
Virucidal activity against enveloped viruses (enveloped viruses including HBV, HCV, HIV as well as coronaviruses)	EN 14476 DWV / RKI guideline	30 sec. 30 sec.
Noroviruses	EN 14476	5 min. ¹⁾

1) Testing under low burden.

Sales

Dürr Dental SE
Höfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
www.duerrdental.com
info@duerrdental.com

Manufacturer

orochemie GmbH + Co. KG
Max-Planck-Straße 27
70806 Kornwestheim
Germany
www.orochemie.de
info@orochemie.com

