AUTOMATIC BLOOD PRESSURE MONITOR

M-PRO

USER MANUAL

dr-frei.com

(EN) **AUTOMATIC BLOOD PRESSURE MONITOR**
M-PRO
User's manual (2-20)

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DEAR DR. FREI AUTOMATIC BLOOD PRESSURE MONITOR OWNER

Thank you for choosing Automatic Blood Pressure Monitor **Dr. Frei**® TM M-PRO. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Trademark **Dr. Frei**®.

Before starting to use this device please study the user manual carefully. The user manual offers all information you need to correctly measure your blood pressure and pulse. For all questions concerning the device please contact your local distributor or **Dr. Frei**® service center in your country.

ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

Type BF applied part.

Read the instructions carefully before using this device.

TABLE OF CONTENTS

What Is Blood Pressure? 4
Which Values Are Normal? 5
Intended Use /Indications for Use 6
Package contents 6
Explanaion of Marks or Symbols 6
Specification 8
Know Your Unit 8
Function setting 9
Tie the cuff 10
Start the Measurement 11
Recall records 11
Delete record 12
Battery maintenance, battery charging 12
Adapter use (optional) 13
Safety Information 13
Tips for operation 15
Common Q&A on blood pressure 15
Care and maintenance 15
Trouble Shooting 16
Electromagnetic Compatibility 16
Periodic safety checks 19
Warranty 20
Disposal 20

IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT


What Is Blood Pressure?

Blood pressure is the pressure that your blood extends to the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a «pump», sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.
There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteriae; and a diastolic (lower) one, which means the blood pressure between two heartbeats.
Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors - time of day, person's condition, physical or mental activity, environment, etc.

Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic. One of the features of the hypertension is its ability to remain unnoticed for

the patient at its early stages. That's why the self-control of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitalises (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

Which Values are Normal?

The world standard as for the norms of the blood pressure is the Classification* of the World Health Organization (WHO):
Table 1. Definitions and classification of blood pressure levels (mmHg)

Category	Systolic	Diastolic
Optimal	< 120	< 80
Normal	< 130	< 85
High-normal	130 - 139	85 - 89
Grade 1 hypertension (mild)	140 - 159	90 - 99
Subgroup: borderline	140 - 149	90 - 94
Grade 2 hypertension (moderate)	160 - 179	100 - 109
Grade 3 hypertension (severe)	≥ 180	≥ 110
Isolated systolic hypertension	≥ 140	< 90
Subgroup: borderline	140 - 149	< 90

When a patient's systolic and diastolic blood pressures fall into different categories, the higher category should apply.

* printed with curtailments

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

If an irregular heartbeat interval was detected in the measurement stored in memory, the icon will be displayed.

ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor.
When measured correctly, if diastolic blood pressure is more than 110 mmHg, it is necessary to call the doctor immediately.

INTENDED USE / INDICATIONS FOR USE

The Arm Blood Pressure Monitor is intended to measure the systolic and diastolic blood pressure as well as the pulse rate of person via non-invasive oscillometric technique in which an inflatable cuff is wrapped around the upper arm. It can be used at medical facilities or at home. The intended upper arm circumference is 22-42cm. Suitable for adults and adolescents who over the age of 12.

PACKAGE CONTENTS

1x Arm Blood Pressure Monitor
1x User Manual
1x Power Adapter
1x Cuff
1x Type-C Cable
1x Travel Bag

EXPLANATION OF MARKS OR SYMBOLS

The below signs might be in the user manual, labeling or other component. They are the requirement of standard and using.

	Follow instructions for use.
	IMPLICATION OF SYMBOL Type-BF applied part.
	Caution: Consult accompanying documents.
	Transport package shall be kept away from rain.
	Do not roll
	Transport package shall not be exposed to sunlight.
	Indicates correct upright position of the transport package.

	Contents of the transport package are fragile therefore it shall be handled with care.
	Indicates temperature limits within which the transport package shall be stored and handled.
SET/	Set User, Voice, Unit
START STOP	Start measurement,Stop measurement, Turn off
MEM/	Memory icon, when the MEM icon flashing, operator could save or recall memories.
IP21	IP21:Protected against solid foreign objects of 12,5 mm ø and greater,Protection against vertically falling water drops.
	Non-ionizing electromagnetic radiation
LOT	Production lot number
	Production date
	Manufacturer
	The device should not be used after the end of the shown or the day
CE 0123	CE mark: indicates that the device complies with the EU 2017/745.
EC REP	Authorized Representative in the European Community
MD	Medical Device


UDI template:
(01):DI(Device Identification)
(17):Expiration Date
(11):Production Date
(10):Lot Number
(21):Serial Number


SPECIFICATION

Name	Automatic Blood Pressure Monitor M-PRO
Model	AES-U522
Measurement mode	Oscillography
Operating conditions	5°C~40°C, 15%-80%RH, 70-106kPa
Storage & transportation condition	-20°C~+55°C, 15% -93% RH, 70-106kPa
Display range	0-290mmHg (0-39kPa)
Measurement range	Diastolic: 30-200mmHg, Systolic: 60-255mmHg Pulse: 40-199 pulses/min
Measurement accuracy	Pressure: ±3mmHg (±0.4kPa), Pulse: ±5% of reading
Memories recall	90 measurement recalls each user
Dimensions	129*109*71mm
Weight	About 270g
Power supply	3,7 V, 800mAh Li-ion battery
Type-c power supply	DC 5V, 1A
Switch off	Automatically turn off after 60 seconds
The cuff size is suitable for arm size	About 22cm between 42cm
Product life	5 years or 10000 measurements under normal use
Adapter	Input:100-240VAC,50/60Hz 0,5A(MAX) output:5V ⎓ 1A which should be applied to IEC60601-1

KNOW YOUR UNIT


A. Number of User Groups
B. Low Battery Symbol
C. High-Pressure Value (Systolic Pressure)
D. Low-Pressure Value (Diastolic Pressure)
E. Blood Pressure Level Indication
F. Cuff Detection
G. Pulse Rate Display
H. Heartbeat Icon


FUNCTION SETTING


Turn off the unit.
Step1: Press the «⚙» button for 3 seconds, when the figure symbol is flashing, means enter the user setting mode.
Step2: Press the «MEM» button, choose user or user
Step3: And then press the «⚙» button and enter the next parameter setting.

Repeat the Step2 and Step3 to set the voice, Set "ON" to open voice and "OFF" to close voice.


Repeat the Step2 and Step3 to set the unit, Set to "ON" in kPa; Set it to "OFF" and the unit is mmHg. The unit will be closed automatically after setting.



TIE THE CUFF

Tie the cuff on your upper arm, and position the tube off-center toward the inner side of the arm in line with the middle finger.

The cuff should not be too tight. You should be able to insert one finger between the cuff and your arm.

These situations could lead to higher blood pressure values due to strain or the arm cuff being lower than the heart.


Sit comfortably with your left arm resting on a flat surface.

Rest for 5 minutes before measuring.
Wait at least 4 ~ 5 minutes between measurements. This allows your blood circulation to recover. For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.


START THE MEASUREMENT

Select the user or user
Press the Start/Stop Button to measure BP automatically.
Press the Start/Stop to turn off, or it will turn off automatically within 1 minute.

RECALL RECORDS

Please choose user or user when the Arm Blood Pressure Monitor turns off, and press the "MEM" to show the average of the last 3 records.

Press the "MEM" to turn over the records.

DELETE RECORD

Delete all memory
In the off state, press the "MEM "button for 3s,
and then press the " ⚙ " button, the display shows " ALL ", which means all the memories have been deleted.


BATTERY MAINTENANCE, BATTERY CHARGING

The blood pressure monitor comes with a rechargeable Li-ion battery inside. Use Type-C USB cable to charge the monitor. will appear on the screen when the monitor needs to be charged.
1. Insert one end of the USB charging cable into the Type-C power socket on the side of the monitor.
2. Insert the other end of the USB charging cable into a DC 5V, 1A adapter outlet, or a powered USB outlet,
3. When the blood pressure level indicator is full grid display, means that the battery is fully charged.

ADAPTER USE (OPTIONAL)

1. The optional AC adapter should comply with the requirements of IEC 60601-1, and all configurations should meet the requirements of the medical electrical system. Please note that local laws take precedence over the above requirements. In case of any doubt, please consult your local representative or technical service department.
2.When using AC power, select the power cord of the TYPE-C jack, and plug it into the power jack.
3.To remove the AC adapter, first unplug the adapter from the outlet and then disconnect the power cord from the device outlet.

Adapter parameter requirements:
Output voltage: DC 5V
Maximum output current: at least 1A
TYPE-C socket
Note: do not insert the adapter when using the product.

⚠ Note

1.Charging may take 3-4 hours.
2.Every 2 months (or when battery life is significantly shorter) , fully charge the monitor and then allow the battery to drain until the monitor shuts off. This will optimize battery performance.
3.Battery life depends on frequency and time of use. If battery life is unusually reduced, please contact Customer Support.
4.Only use a Type-c USB charging cable tocharge the monitor.
5.In extreme conditions, the battery may leak corrosive fluid. If this comes into contact with eyes or skin, rinse immediately with water and seek medical attention. Do not replace the lithium battery by yourself, inadequately trained personnel could result in a hazard.

SAFETY INFORMATION

1. General usage
Do not adjust medication based on measurement values from this blood pressure monitor. Take medication as prescribed by your physician. Only a physician is qualified to diagnose and treat High Blood Pressure.
- The monitor is not intended to be a diagnostic device.
- Consult your physician before using the device for any of the following conditions: common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, preeclampsia, and renal diseases.
- Note that PATIENT motion, trembling, and shivering may affect the measurement.
- Do not use the device on the injured arm or the arm under medical treatment, as this

can cause further injury.
- Do not apply the arm cuff on the arm while on an intravenous drip or blood transfusion.
- Prolonged over-inflation of the monitor will result in harmful injury to the patient.
- Too frequent measurements can cause injury due to blood flow interference.
- Consult your physician before using the device on the arm with an arterio-venous (A-V) shunt.
- Do not use the device with other medical electrical (ME) equipment simultaneously.
- Do not use the device in the area of HF surgical equipment, MRI, or CT scanner, or in an oxygen-rich environment.
- Please ask your doctor about your normal blood pressure for the right direction before taking the measurement by yourself.
- If the cuff causes any discomfort, please turn off the equipment by pressing the START/STOP button.
- If the arm cuff doesn't deflate automatically after the equipment has pressurized to 300mmHg (40kPa), please press the START/STOP button to stop inflation and take off the cuff.
- This product applies only to adults. Please keep the unit out of reach of children.
- The device complies with RF specifications when the device is used at 0mm from your body.
- This monitor is calibrated at the time of manufacture, if the monitor is used according to the instruction, periodic recalibration is not required. If it is inaccuracy often, please contact your retailer or customer services.
- Do not disassemble, repair, or remodel the main unit or the cuff of the blood pressure monitor by yourself. If necessary, contact your retailer or customer services.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
- The patient can use and maintain the device as an operator.
- Please do not continue to use the cuff if you are allergic to it.
- Please use our company's matching cuff, if you use other parts or materials may not be able to boot or reduce safety.
- If non-manufacturer-supplied parts are used, errors in measurement results may be caused.
- Please check whether the battery is leaking before use, if it leaks, do not use it.
- Please check whether the power adapter is in good condition before use. If it is damaged, please replace it with a new one.
- Do not remove the power adapter during the measurement when using the power adapter to supply power.
- The operator shall not touch the output of the AC adapter and the patient simultaneously.

2. General precautions
- Do not forcibly crease the arm cuff or the air tube excessively.
- Do not press the air tube while taking a measurement.
- Do not drop the monitor or subject device to strong shocks or vibrations.
- Do not inflate the arm cuff when it is not wrapped around your arm.
- Do not use the device outside the specified environment. It may cause an inaccurate reading.

- Dispose of the device, components, and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.

TIPS FOR OPERATION

It can cause incorrectness if the measurements are taken in the following circumstances.
- Less than 1 hour after eating or drinking.
- Immediate measurement after tea, coffee, and smoking.
- Less than 20 minutes after taking a bath.
- When talking or moving your fingers.
- In a very cold environment.
- When you want to discharge urine.

COMMON Q&A ON BLOOD PRESSURE

Q1: Why does my blood pressure vary even in one day?
1. Individual blood pressure varies even in one day, it is also affected by the way you tie your cuff and your measurement position, so please take the measurement in equal conditions.
2.The difference of the pressure is greater if the person takes medicine.
3.Wait at least 4 ~ 5 minutes between measurements.
Q2: Why the blood pressure I get from the hospital is different from home?
The blood pressure is different even during 24 hours because of the weather, emotion, exercise, etc.
Attention is required.
When you measure your blood pressure at home:
If the cuff is tied properly.
If the cuff is too tight or too loose.
If the cuff is tied on the upper arm.
If you feel anxious or pressured.
Take deep breaths 2-3 times before beginning.
Advice: adjust yourself for 4-5 minutes until you calm down.
Q3: If the result is the same if measuring on the right arm:
It is ok for both arms, but there will be some different results for different people, so we suggest you measure the same arm every time.

CARE AND MAINTENANCE

1. Cleaning the monitor and cuff
- Make sure the monitor is off prior to cleaning, it can be sterilized with 70%-80% medical alcohol on the soft towel or cotton sliver.

- Clean the monitor with a soft dry cloth, Do not use any abrasive or volatile cleaners.
- Use a soft moistened cloth and soap to clean the arm cuff, do not wash it in water.
- Never immerse the monitor or any of the components in water.

2. Maintaining the monitor and cuff
- Keep the monitor and cuff in the storage box when not in use.
- Do not forcefully bend the arm cuff or air tube. Do not fold tightly.
- Do not press the START/STOP button before wearing the cuff properly.
- Do not disassemble or attempt to refit the unit or components.
- Do not subject the monitor to strong shocks, such as dropping the unit on the floor.
- Protect the unit from contamination, dust, and direct sunlight.
- If the unit will not be used for a long time, please remove all the batteries.

TROUBLE SHOOTING

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	The display is dim or will not light up	Batteries are exhausted	Recharge the battery
Low batteries		Batteries are low	Recharge the battery
Error message	Err 1	The heart beat signal is too weak or the pressure of cuff decreased	Refasten the cuff and thenmeasure again.
	Err 2	There are some noise to affect measurement.	Refasten the cuff and then measure again.
	Err 3	The measurement is not right.	Refasten the cuff and then measure again.
	Err P	Inflate failed	Refasten the cuff and then measure again
	Err H	Overcharging pressure	Refasten the cuff and then measure again

ELECTROMAGNETIC COMPATIBILITY

The ME EQUIPMENT or ME SYSTEM is suitable for a home or hospital environment.
Warning: Don't use near active HF surgical equipment and the RF-shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.
Warning: Use of this equipment adjacent to or stacked with other equipment should

be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
If any: a list of all cables and maximum lengths of cables (if applicable),
transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g. shielded cable, load impedance) or specifically (e.g. by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).
If any: the performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term "ESSENTIAL PERFORMANCE" need not be used).
1, All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
2, Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacturer´s declaration – electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations /flicker emissions IEC 61000-3-3	Applied

Table 2

Guidance and manufacturer´s declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Power supply lines: ±2 kV
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV, line(s) to earth: ±2 kV. 100 kHz repetition frequency	line(s) to line(s): ±1 kV. 100 kHz repetition frequency

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0,5 cycle At 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 ° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	0% 0,5 cycle At 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 ° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

NOTE U T is the a.c. mains voltage prior to application of the test level.

Table 3
Guidance and manufacturer's declaration - electromagnetic Immunity

WARNING
Do not modify this equipment without authorization of the manufacturer.
All major maintains on the device must be performed by an authorized service center or distributor. No use-serviceable parts inside, before servicing to authorized representative or manufacturer!

DECLARATION:
Do not modify this equipment without authorization of the manufacturer.
All major maintains on the device must be performed by an authorized service center or distributor. No use-serviceable parts inside, before servicing to authorized representative or manufacturer!

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSUR PORT IMMUNITYto RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27
	450	380 – 390	GMRS 460, FRS 460	FM ±5kHz deviation 1 kHz sine	2	0.3	28
	710 745 780	704– 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	810 870 930	800– 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	1720 1845 1970	1700– 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3,4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	2450	2400– 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240 5500 5785	5100– 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9

PERIODIC SAFETY CHECKS

If you use the device with power adapter, preventive inspection and maintenance to be performed including the frequency of such maintenance.
- Every time before use, please check the adapter, once damaged, never to use.
- Please clean the plug of adapter plug at least once a year. Too much dust on plug may cause the fire.

The manufacturer reserves the right to make technical changes without notice in the interest of progress.
Prior notices will not be given in case of any amendments within this manual. The mentioned trademarks and names are owned by the corresponding companies.

WARRANTY

Your Automatic Blood Pressure Monitor M-PRO is warranted for 5 years from date of purchase. Warranty for the cuff is 1 year from the date of purchase. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.
The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed.

Disposal
Do not dispose of electrical appliances as unsorted municipal waste, use separate collection facilities. Contact you local government for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

Manufacturer information:
Manufacturer: Alicn Medical Shenzhen, Inc
Address: Room 4C, Building 1, Zhongcheng Biomedical Industrial Park, No. 21 Linhui Road, Jinsha Community, Kengzi Street, Pingshan District, 518118 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA
Tel: 0086-755-26501548 Fax: 0086-755-26504849
Email: Info@alicn-med.com Website: http://www.alicn-med.com

EC REP **Authorized EU Representative**
Company name: CMC Medical Devices & Drugs S.L.
Address: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain
Web: www.cmcmedicaldevices.com
Tel: +34951214054 Email: Info@cmcmedicaldevices.com

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