AUTOMATIC BLOOD PRESSURE MONITOR

M-MAX

USER MANUAL

dr-frei.com

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DEAR DR. FREI AUTOMATIC BLOOD PRESSURE MONITOR OWNER

Thank you for choosing Automatic Blood Pressure Monitor **Dr. Frei®** TM M-PRIME We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Trademark **Dr. Frei®**.
Before starting to use this device please study the user manual carefully. The user manual offers all information you need to correctly measure your blood pressure and pulse. For all questions concerning the device please contact your local distributor or **Dr. Frei®** service center in your country.

ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

Type BF applied part.


Read the instructions carefully before using this device.

EN

TABLE OF CONTENT

Parts and components 6
Symbols 6
General 7
Important safety instructions 8
Use power adapter 12
Setting date and time 12
Voice function adjustment 13
Setting the languages 13
Use the device 14
Correct posture 14
Assembly the cuff 15
Carry out a measurement 16
Function of memory 17
Who classification indication 18
Error and low battery information 19
Troubleshooting 20
Care, storing, repairand recycling 20
Specification 21
Manufacturer's declaration 22
Periodic safety checks 25
Warranty 25

PARTS AND COMPONENTS

1. Main Body
2. Display
3. Air Connector
4. Tube Plug
5. Air Hose
6. Cuff
7. D-ring
8. Button ‘ ⏻ ’
9. Button ‘ M1/M2 ’
10. Adapter
11. Type-C cable
12. Batteries
13. Travel Bag

SYMBOLS

SYMBOLS	Symbols Meaning
	Manufacturer
	Authorized Representative in the European community

	Symbol for the marking of electrical and electronics devices according to Directive 2012/19/EU. The device, accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal.
CE 0123	CE marking in conformity with EU directive 2017/745
	Keep dry
	Attention, consult accompanying documents
	Type BF Applied Part
	Stand by
	CLASS II
MD	Medical Device

GENERAL

This instruction manual is intended to assist the user for safe and efficient operation of the automatic digital talking blood pressure monitor (hereinafter: device) model M-MAX. The device must be used in accordance with the procedures described in the manual. It is important to read and understand the entire manual, especially the section < IMPORTANT SAFETY INSTRUCTIONS >. This device is intended for the non-invasive measurement of systolic and diastolic arterial blood pressure and pulse rate in adults (age 15 and above).

CAUTION:
1. Do not use this device on infants or persons who cannot express-their intentions.
2. The device is not suitable for measuring the blood pressure of children.

Ask your doctor before using it on elder children.
3. The patient is an intended operator. But persons who suffer from arrhythmia, diabetes, cardiovascular problems or who have had a stroke should consult your doctor before using the device.

PRINCIPLE OF OPERATION

This device adopts the oscillometric technology with Fuzzy Algorithm to measure the arterial blood pressure and pulse rate. The cuff is wrapped around the arm and automatically inflated by the air pump. The sensor of the device catches weak fluctuation of the pressure in the cuff produced by extension and contraction of the artery of the arm in response to each heartbeat. The amplitude of the pressure waves is measured, converted in millimeters of the mercury column, and is displayed by digital value.

ATTENTION: This device can not provide reasonable accuracy if used or stored in the temperature, humidity or altitude beyond the range stated in the section <SPECIFICATIONS> of this manual.

NEW TECHNOLOGIES USED

Fuzzy Algorithm is the processing algorithm, taking into account the specialty of Individual heartbeats; which provides higher accuracy of measurement." Software version: V1.1

IMPORTANT SAFETY INSTRUCTIONS

It is necessary to know that arterial blood pressure is subjected to sharp fluctuations.The level of the arterial blood pressure depends on many factors. Generally arterial blood pressure is lower in summer and higher in winter. Arterial blood pressure changes with atmospheric pressure and is affected considerably by many factors, e.g. physical loads, emotional excitability, stress, meals, etc. Medicines, drinking, smoking affect greatly the level of an individual's blood pressure. Blood pressure does vary with age and individuals, and it is recommended to write down the readings from blood pressure records daily, then you can check with your doctor to find out what is a "normal blood pressure measurement" for you.

Please read the instruction manual carefully before using this device, especially < Important safety instructions^ it can help you use the device correctly and safely!

Please keep the instruction manual for future use. For specific information about your own blood pressure, consult your physician.

Warnings

• When pets or children play, the air tube may be bent or the device may be accidentally dropped, affecting the product performance. Therefore, please protect the device to prevent pets and children from affecting the

use of the product when use.
• Consult your physician if you suffer from illnesses prior to using the device.
• The device is not suitable for persons who have electrical implants.
• If you had a mastectomy (breast amputation) do not use this blood pressure monitor on the arm on the side of the mastectomy.
• Pregnant women should only measure their own blood pressure in consultation with their doctor, since the readings may be changed with pregnancy.
• Do not service or maintain the cuff while in use with patient.
• Do not use this blood pressure monitor on any arm where intravascular access or therapy (such as an intravenous drip or a blood transfusion), or an arteriovenous shunt (A-V shunt) is present. The temporary interference to blood flow by the blood pressure measurement could result in injury.
• Do not use the device with other medical electrical (ME) equipment simultaneously.
• Do not use the device in the area the HF surgical equipment, MRI, or CT scanner exists, or in the oxygen rich environment.
• Do not use a mobile phone or other devices that emit electromagnetic fields, near the device. This may result in incorrect operation of the device.
• Never use any accessories or parts from other manufacturers. Using such accessories or parts could cause a hazardous situation for the user or damage to the device.
• Not modify this equjpment without authorization of the manufacturer
• The batteries used in this device may present a fire or chemical burn hazard if mistreated. Do not disassemble, heat or incinerate.
• Keep equipment away from fire and heat sources to prevent fire or explosion
• Please keep the unit out of reach of infants, children or pets, since inhalation or swallowing of small parts can be dangerous or even fatal.
• Please pay attention that the continuous CUFF pressure due to connection tubing kinking will cause a harmful injury.
• Do not use an extension cord with this device.
• The air tube or the AC adapter cable may cause accidental strangulation in infants.
• Do not put the air tube around your neck - danger of suffocation!
• The device should never be left unattended when plugged in.
• Do not reach for a corded device that has fallen into water. Unplug immediately.
• It is quite normal that two measurements taken in quick succession may produce significantly different results, because too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.

Cautions

• Use this device under the right environmental conditions as indicated in

this user manual. If not, this could affect the performance, lifetime of the device and measurement results.
• Only use this device for its intended purpose as described in this user manual.
• Do not confuse self-monitoring with self-diagnosis. This device allows you to monitor your blood pressure. Do not begin or end medical treatment based on the
measurement results. Always consult your physician for treatment advice.
• Do not take any therapeutic measures on the basis of a self-measurement. Never change prescribed medication without consulting your physician. Consult your physician if you have any questions about your blood pressure.
• If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure.
• Consult the physician if measurement errors occur in children or persons with arrhythmia.
• The pulse display is not suitable for monitoring the frequency of cardiac pacemakers.
• Common arrhythmias (such as atrial or ventricular premature beats or atrial fibrillation) and peripheral artery disease / arteriosclerosis can affect the accuracy of this blood pressure monitor. Please consult your physician how to best use this blood pressure monitor if you suffer from any of these conditions. Blood pressure measurement is not suitable in cases of serious arteriosclerosis (hardening of the arteries).
• The effectiveness of this blood pressure monitor has not been established in pregnant women.
• Always check the device and cuff before you use it. Do not use the device or cuff if one of them is damaged, because this may cause injury.
• This device is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
• Do not attach the cuff on the same arm on which other monitoring medical electrical equipment is attached simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring medical electrical equipment.
• Never attach the cuff on injured skin, an injured arm or an arm under medical treatment as this can cause further injury.
• Do not forcibly crease the arm cuff or the air tube excessively.
• Do not press the air tube while taking a measurement.
• Do not use the device in case of existing polyester or nylon material allergies.
• This device is not suitable for continuous monitoring during medical emergencies or operations.

• This device cannot be used with HF (High Frequency) surgical equipment at the same time.
• This device is not washable. Never immerse the device in water and do not rinse it under the tap.
• This device should keep dry to prevent from moisture.
• The equipment is not AP/APG equipment and is not suitable for use in the presence of a flammable anesthetic mixture with air, with oxygen or nitrous.
• To avoid measurement errors, do not use the device near strong electromagnetic fields, radiated interference signal or electrical fast transient/burst signal. For example magnets, radio transmitters, microwave ovens.
• If this device was stored in low temperature, leave it in room temperature for at least 1 hour.
• Repeated measurements with an interval of 3 minutes are recommended, so you can calculate the average to get a more accurate measurement. An internal of 3 minutes can also ensure that the operation of the device does not result in prolonged impairment of the circulation of the blood.
• Atherosclerosis patients may require longer interval (10-15minutes) as elasticity of patient's vessels decreases significantly with the disease. 10-15minutes interval is also applicable for patients suffering from diabetes for a long period of time.
• Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
• Connecting electrical equipment to mso effectively leads to creating a ME system, and can result in a reduced level of safety.
• Ensure that this monitor has accliamed to room temperature befor measurement. Taking a measurement after an extreme temperature change ! could lead to an inaccurate reading.We recommends waiting for approximately 2 hours for the monitor to warm up or cool down when the monitor is used in an environment within the temperature specified as operating conditions after it is stored either at the maximum or at the minimum storage temperature.

CLASSIFICATION

• ME EQUIPMENT not intended for use in an oxygen rich environment or in the presence of flammable mixers.
• Internally powered equipment (without adaptor), Class II equipment (with adapter).
• Type BF applied part, recognize the cuff as applied part.

USE POWER ADAPTER


Besides batteries you can use power adapter as the power supply. The power adapter is optional for this device for sale.
The adapter is specified as a part of the blood pressure monitor.
• Insert the adapter cord into the jack on the right side of the monitor
• Insert the adapter plug into the outlet.
• To remove the adapter, disconnect the adapter plug from the outlet first and then disconnect the cord from the monitor's jack.

Caution

• When using optional adapter, the adapter must comply with the requirements of standard IEC60601-1.
• To avoid possible damage of the monitor, use only the exclusive adapter that can be purchased from authorized dealers. Other adapter may damage the blood pressure monitor.
• The adapter is used as an isolating means, the adapter plug shall insert into the outlet nearby the operator, make it easy to disconnection the device from the outlet.
• If long time work, remove the plug after the adapter cools, and prevent bums.
• Plug the adapter into the appropriate voltage outlet. Do not use in a multi-outlet plug.
• Not to position the blood pressure monitor to make it difficult to operate the disconnectionLdevice(adaptor).

Note: The monitor is designed not to draw power from the batteries when the adapter in use.

Adapter technical feature:
Model: YS5M-0501000
Input: 100-240V 50/60Hz
Output voltage: 5V
Output current: 1A

SETTING DATE AND TIME

The function provides accurate measuring time for each measurement. To get an accurate date and time, the user should preset the date and time correctly before the first use of this device. The operation procedure for presetting Date/Time is as follows:
1. When the device is on stand by, press and hold the button 'M1', then press the button ' ⏻ ' , and the year number flashes as in Fig. 1;

Fig.1

2. Press button ‘M1’ to add the number, and press button ‘M2’ to decrease the number,then press the button ‘ ⏻ ’ for confirmation;
3. When the year setup is finished, the month number will flash automatically. Please follow the same instruction as above to set month, date and time as in Fig.2;
4. If you want to change the date and time, please repeat procedure 2,3,4.
Annotation: When under the mode of Date/Time setting without any operation, the device will automatically save the operation and return to shutdown mode after 30s.

Fig.2

VOICE FUNCTION ADJUSTMENT

The device has a voice function. To adjust the on/off of the voice, please follow procedures below:
1. When the device is on stand by, press and hold button “M1” , then press button “ ⏻ ” two times until “ON” or “OF” will display on the LED and “ 🔇 ” flash as in Fig.3;
2. Press button ‘M1 to adjust the voice on or off;
3. Press button “ ⏻ ” for confirmation.
Annotation:
1. When the device is used for the first time,the voice function is enabled by default.
2. When under the mode of VOICE FUNCTION ADJUSTMENT without any operation, the device will automatically save the operation and return to shutdown mode after 30s.

Fig.3

SETTING THE LANGUAGES

There are five languages of voice for choice. To adjust the languages:
1. When the device is on stand by,press button ‘M1’ and don’t release, then press button ‘ ⏻ ’ three times until the language icon displays on the LED as in Fig.4;
2. The device will show “L1” (English is the default when the device is used in first time)with the voice prompt. Press button ‘M1’ to set the other languages ;
3. Press button ‘ ⏻ ’when you have finished setup. If you want to change the languages, please repeat the above procedure. L1-English L2-Bulgarian L3—Romanian L4—Russian L5-Spain

Fig.4

Annotation:
1 .The icon “L1” in the LED indicates the language is English, And there will be

a demo voice “We wish you good health.” in English.Other icons “L2-L5” speak a demo voice in relative languages.
2. The voice function will increases the power consumption and reduces the life of batteries.
3. When under the mode of SETTING THE LANGUAGES without any operation, the device will automatically save the operation and return to shutdown mode after 30s.


USE THE DEVICE

Caution:
• Please keep quiet for 5-10 minutes and avoid eating, drinking, alcohol, smoking, exercising and bathing before taking a measurement. All these factors will influence the measurement result.
• Remove any garment that fits closely to your upper arm.
• Always measure on the same arm (normally left).
• Measurements should be taken regularly at the same time of each day, as the blood pressure varies even during the day.
• Any effort to support the arm during measurement may increase the measured blood pressure.
• Make sure, you are in a comfortable, relaxed position with leg-uncrossed, feet flat on the floor, back and arm supported, middle of the cuff at the level of the right atrium of the heart and do not move or constrict your muscles and talk during measurement. Use a cushion to support your arm if necessary. Keep position in normal use.
• If the arm artery lies lower or higher than the heart, a false reading will be obtained.
• A loose or open cuff causes false readings.
• With repeated measurements, blood accumulates in the arm which can lead to false reading.
• Consecutive blood pressure measurements should be repeated after 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away.

CORRECT POSTURE

1. Sit beside the table and let the table support your arm as you take the measurement.
2. Sit upright with your back straight.
3. Make sure that the cuff on the upper arm no cross, and is at approximately the same level as the heart.
4. Make sure that your feet lie on the ground and no cross.
5. You may lie on your back and take a measurement. Look at the ceiling, keep calm, and don’t move your neck or body during the measurement.

ASSEMBLY THE CUFF

1. Insert the edge of the cuff approximately 5 centimeters Into the D-rlng as shown in figure.
2. Put the cuff on the left upper arm with the tube pointing to the direction of palm. If measurement on your left arm is difficult, you can use right arm fbr measurement. In this case, it is necessary to know that the readings may differ about 5-10 mmHg between left arm and right arm.
3. Wrap cuff around your upper arm with the lower edge of the cuff approximately 2-3 centimeters above the elbow. The mark <ARTERY> must be over the artery of the arm.
4. Press the cuff to make sure that it is attached securely. The cuff should not be too tight or loose is greatly recommended. Two fingers should be easily put in between cuff and upper arm.
5. The mark <INDEX> on the cuff must point to area <NORMAL>. This means the cuff size is correct. If mark <INDEX> points to the area beyond area <NORMAL>, please consult your dealer whether you need another
size cuff. This device is supplied with the standard cuff which is fit for the arm size 22-42 cm.
6. Sometimes it is difficult to make the cuff regular depending on the shape of the user's upper arm, the cone-shape assembly of cuff is also acceptable.
7. If your clothes restrict the blood circulation of your upper arm, or you roll your sleeve up so as to result in such restriction. Please take off your shirt to get an accurate measurement if necessary.

Caution:

If you experience discomfort during a measurement, such as pain in the upper arm or other complaints, press the ' ⏻ 'button to release the air immediately from the cuff. Loosen the cuff and remove it from your arm.

CARRYING OUT A MEASUREMENT

1. Insert the air hose into the air connector. Before the measurement, take 3-5 times deep breath and relax yourself. Don't talk or move your arm.In off state,press button ' ⏻ ' to enter the standby, the LED displays the time and date alternately.

Fig.5

2. When the device is on stand by,press button ' ⏻ ', and all symbolswill appear on display as Fig.5 with voice prompt: "Please relax and don't move, we'll start to measure now.". (Press button "M" to skip the voice prompt) Then '0 mmHg' will appear on the screen. Pump begins to inflate with display showing the reading of pressure.

Fig.6

3. If the cuff was applied too loosely, it may cause unreliable measurement results or measurements can fail to start. The " **Cuff Check Indicator** " will flash to determine if the cuff is wrapped snugly enough. The specified icon ' ' appears if the cuff is wrapped correctly during measurement as Fig.6. Otherwise machine error appears if the cuff is wrapped too loosely see the error information for details.

Fig.7

4. Pressure in cuff increases to working pressure. Then the pump stops and pressure falls quickly.

NOTE: The device will inflate to a higher pressure automatically in case the inflation pressure is not sufficient to determine measurement result.

5. The heart symbol " " start to flash every heart beat once pulse is detected as Fig.7.

Fig.8

6. After the measurement, the ' ' or ' ' will flash to remind the user to record the reading, and systolic pressure, diastolic pressure, pulse rate as Fig.8. Press button 'M1 or 'M2' to record the reading in corresponding memory.The device will shut down and the result will not be saved if no operation is performed within 30 s.

Fig.9

7. If irregular heartbeat was detected during the measurement, LED display the ' ' icon to remind users of heartbeat irregularity as Fig.9.

Attention: We recommend contacting your physician if you see the ' ' indicator frequently.

8. The ***"Movement Detection"*** helps reminding the user to remain still and is indicating any adverse body movement during measurement. The specified icon ' ' appears once a "body movement" has been detected

during and after such a measurement.
Note: It's highly recommended that you measure again if the icon appears.
9. Press the button ' ⏻ ' to turn off the device. Please rest for at least 3 minutes for another measurement.

NOTES:

If the result exceeds the display range,' Hi ' is displayed. When the result is below the display range, ' Lo ' is displayed.
RAPID DEFLATION DURING MEASUREMENT
If you do not feel well during measurement or want to stop the measurement for some reason, you can press the Button " ⏻ " .The device will quickly release the air in cuff and the device will be switched off.

FUNCTION OF MEMORY

MEMORY RECALL

1.The Device can store 100 sets of readings each in 'M1' and 'M2', and automatically calculate the average value of the latest 3 readings. When the memory is full (100 sets of readings are stored), the oldest reading will be replaced by new one automatically. Memory will not clear away even if power supply is removed; 2.After a measurement or when the device stands by, the user can press button 'M1' or 'M2' to recall memory. Press button 'M1°, the display will show the average value of the latest 3 readings with a voice prompt as Fig. 10; 3.Press again, the display will show '01' which means the latest reading and measuring time with a voice promptas,then turns to another screen to show readings Fig. 11; 4. Press again, the display will show '02', which means the second to the llatest reading...

Fig. 10


Fig. 11

MEMORY CLEARANCE


When the device is switched off or when the device stands by, hold down button 'M1' for at least 3 seconds, the display will show 'CL' which means all the stored reading are removed as Fig.12. During the memory query,hold down button 'M1' or 'M2' for at least 3 seconds, the display will show 'CL' which means all the stored reading in this group are removed as Fig.13.

Fig. 12 Fig. 13

WHO CLASSIFICATION INDICATION

Standards for assessment of high or low blood pressure, regardless of age, have been established by Work! Health Organization (WHO) as show in the chart as below:

Arterial blood pressure classification by World Health Organization

The indicator display a segment, based on the current data, corresponding to the
WHO classification.
For example, if your blood press is 145mmHg (systolic pressure), 88mmHg (diastolic pressure), according to the world health organization standard, your blood pressure level is Mild Hypertension.
Note: If the systolic blood pressure and the diastolic blood pressure fall into different categories, the higher value should be taken for classification.

IRREGULAR HEARTBEAT DETECTOR

Model M-MAX digital blood pressure monitor provides a blood pressure and pulse rate measurement even when an irregular heartbeat occurs. When the device detects the irregular heartbeat or any excessive body movement during measurement, the ' ' icon will display in the LCD. It is important that you are relaxed, remain still and do not talk during measurement.
Notice: We recommend contacting your physician if you see this '' indicator frequently.

ERROR AND LOW BATTERY INFORMATION

INDICATION	POSSIBLE REASON	CORRECTION METHODS
	ER00: The device can't find zero point.	Contact the seller or manufacturer
	ER01:Blood pressure and pulse cannot be measured. ER02:The cuff is not wrapped snugly enough. ER03:There is an adverse body movement.	Press button ' ⏻ ' to return the standby.Make sure that cuff is put on correctly and the tube plug is inserted tightly and repeat the measurement. Repeat the measurement with following completely recommendations of manual.
	The batteries are low, the device can still operate normally.	Prepare for battery replacement
	The batteries are weak.	Replace all batteries.

TROUBLESHOOTING

SYMPTOM	CHECK	POINT REMEDY
No display when connect the power.	The batteries have run down.	Replace all batteries
Inflation stops and reinflate later.	The automatic inflation for ensuring correct measurement. Did you talk or move your arm (or hand) during measurement?	Keep quiet and silent during the measurement.
The reading is extremely low or high.	Is the cuff at the same level as the heart? Is the cuff wrapped right? Did you strain your arm during measurement? Did you talk or move your arm (or hand) during measurement?	Make sure that your posture is right. Wrap the cuff correctly. Relax during measurement. Keep quiet and silent during the measurement.
Pulse rate is too low or too high.	Did you talk or move your arm (or hand) during measurement? Did you make measurement right after exercise?	Keep quiet and silent during the measurement. Take measurement again after resting for more than 5 minutes.

CARE,STORING, REPAIRAND RECYCLING

1. It's necessary to protect this device against high moisture, direct sunlight, shock, solvent, alcohol and gasoline.

2. Remove the batteries if the device is being stored for along time, and keep the batteries far away from children.
3. Keep the cuff away from sharp objects and don't extend or twist the cuff.
4. This device is not washable. Never immerse the device in water and do not rinse it under the tap. Use only soft and dry cloth to clean the device.
5. Do not serve or maintain the cuff and the device when in use with patient.
6. The cuff is sensitive and must be handled with care. You can clean the cuff with damp cloth for daily maintenance. To avoid cross infection when sharing the cuff, you can sterilize the fabric cover of the cuff with tampons moistened by 3% solution of hydrogen dioxide. After long use there will be a partial discoloration on the fabric surface of the cuff. Do not laundry the cuff as well as ironing with a hot flatiron.
WARNING: Under no circumstances may you wash the inner bladder!
7. Since neither the device nor batteries are household waste, follow your local recycling rules and dispose them at an appropriate collection site.
8. Do not open the device, or delicate electrical components as an intricate air unit could be damaged. If you can not fix the problem using the troubleshooting instruction, please request service from your dealer.
WARNING: Do not repair the device without manufacturer's authorization. Do not carry out maintenance when using the device.
Caution:
Generally, we recommend the device should be inspected every 2 years and utilize the manometer mode to verify the accuracy of the manometer at least at 50mmHg and 200mmHg after maintenance and repair. Please contact your dealer for maintenance.

SPECIFICATIONS

Model	M-MAX (LD-587)
Size	106(L) *140 (W)x48(H)mm
Weight	Approximately 240g
Measuring method	Oscillometry
Measuring range	40 to 180mmHg(DIA,diastolic pressure) 60 to 260mmHg(SYS,systolic pressure) 40 to 160 beats/minute (PUL,pulse rate)

Measuring accuracy	± 3 mmHg for static pressure ± 5% of the reading for the pulse rate
Inflation	Automatic by the pump
Rapid defiation	Automatic electronic valve
Power	Type-C adapter DC5v, 1A Battery 4XAAA
Memory	2 Users with 100 sets of memory each
Operation temperature and humidity, air pressure	+10°C to + 40 °C, 85% and below 800hPa to 1060hPa
Transport and storage temperature and humidity, air pressure	-20 °C to + 50 °C, 85% and below 500hPa to 1060hPa
Upper arm circumference	Applicable for arm circumference 22-42cm
Complete kit	Main body, cuff, 4*AAA batteries, adapter, cable Type-C, gift box, instruction manual, warranty card, travel bag.
Overvoltage category	Category II

HONSUN (NANTONG) Co., Ltd.
Address:No.8 Tongxing Road Economic & Technological De velopment

EU REP SHANGHAI INTERNATIONAL HOLDING CORP. GMBH (EUROPE)
Address: Eiffestrasse 80, 20537 Hamburg, Germany

MANUFACTURER'S DECLARATION

Compliance information for each EMC test

Electromagnetic Emission(Home Healthcare Environment)	
Emission test (IEC60601-1-2:2014)	**Compliance**
Conducted and radiated RF emissions	CLSPR 11 Group 1 Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies

Compliance information for each EMC test

Declaration-Electromagnetic Immunity(Home Healthcare Environment)		
Immunity test	**IEC 60601 test level**	**Compliance level**
Conducted RF IEC 61000-4-6	3V 150 kHz to 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz	3V 150 kHz to 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz also meet the requirement of table 9 of 60601 -1-2:2014	10 V/m 80 MHz to 2.7 GHz also meet the requirement oftable 9 of 60601-1-2:2014
Electrostatic discharge (ESD) I EC 61000-4-2	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air

Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines
Surge IEC 61000-4-5	±0.5 kV,± 1 kV line(s) to lines	±0.5 kV,± 1 kV line(s) to lines
Voltage dips short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% Ut , 0.5 Cycle at 0°,45°,90°, 135°, 180°, 225°, 270°, 315° 0% Ut , 1 Cycle and 70% Ut, 25/30 cycles sigle phase:at 0° 0% Ut ,250/300 cycles	0% Ut , 0.5 Cycle at 0°,45°,900, 135°, 180°, 225°, 270°, 315° 0% Ut , 1 Cycle and 70% Ut , 25 cycles sigle phase:at 0° 0% UT,250cycles
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m

NOTE: The EUT is the a.c. mains voltage prior to application of the test level. The following phenomenon is still fulfill the requirement of basic safety and essential performance.
*UT:230V ~/50Hz.The pressure of the EUT is deviation the normal value but the value is still more than 10psi when flow is 4.5l/min.
**UT:230V ~/50Hz.The EUT stop working when adding 0%UT,but the EUT can restore its normal mode automatically.
• Use of this equipment adjacent to or stacked with other equipment should be avioded because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
• Portable RF communications equipment(including peripherals such as antenna cables and external antennas) should be used no closer than 30cm(12 inches to any part of this devie, inluding cables specificed by the manufacturer.
Otherwise, degradation of the performance of this equipment could result.
• Under the test condition specified in immunity, the product can provide the basic safety and essential performance.
• If the essential performance is lost or degraded,additional measures are necessary,such as reorienting or relocating the device.
1. Warranty for this automatic digital blood pressure monitor is 24 months from the date of purchase. The 24 months warranty excludes the monitor cuff. The cuff is warranted for 12 months.
2. The warranty obligations are prescribed for by warranty certificate for buyer.

3. The addresses of organizations for guarantee maintenance are present in the warranty certificate.

WARNING

Do not modify this equipment without authorization of the manufacturer. All major maintains on the device must be performed by an authorized service center or distributor. No use-serviceable parts inside, before servicing to authorized representative or manufacturer!

DECLARATION:

When technical information for user or service personnel requirements is not in the scope of confidentiality of the Company, the Company committed to provide information disclosure in accordance with procedure, including circuit diagrams and parts lists, and other related type technology information that do not involve commercial secrets may be disclosed. Access to information channels and procedures, please contact your dealer or manufacturer.

PERIODIC SAFETY CHECKS

If you use the device with power adapter, preventive inspection and maintenance to be performed including the frequency of such maintenance.

• Every time before use, please check the adapter, once damaged, never to use.

• Please clean the plug of adapter plug at least once a year. Too much dust on plug may cause the fire.


The manufacturer reserves the right to make technical changes without notice in the interest of progress.

Prior notices will not be given in case of any amendments within this manual. The mentioned trademarks and names are owned by the corresponding companies.

WARRANTY

Your Automatic Blood Pressure Monitor M-MAX is warranted for 5 years from date of purchase. Warranty for the cuff is 1 year from the date of purchase. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.

The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed.

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