

Tendo®

# must™

multi use single technique





# Shape with Single Use Multiple

You can complete your treatments effectively in a short time with T-endo must, which is designed to eliminate the problems we may encounter in root canal shaping. T-endo must is designed for multiple use.

Therefore, it allows you to save significantly on your treatment costs.

Thanks to the powerful files produced with TM-Wire Heat Treatment technology, a short and effortless preparation will be sufficient for the treatment. You will be amazed at the result when the taper angle designed to protect pericervical dentin is combined with the power of the Reciprocal movement.

Sterile blister packaging offers an ideal solution for the safety of both you and your patient. This ready-to-use package does not require sterilization when first opened.

The 'S' section tips found in T-endo must files are unnecessary in the root canal. Prevents substance loss and provides a more controlled treatment opportunity

"tg Glide Path File", designed especially for narrow channels, A fast and effective reciprocal experience thanks to its square-section structure presents. With these features, T-endo must, root canal shaping It offers a reliable and effective solution.



Sterile



Heat  
Treatment




Reciprocal  
System



Fast




## Cyclic Fatigue Test



Compatible with the world's leading reciprocal file brands.  
NTUB, Number of turns until broken

## Torsional Resistance Test



Compatible with the world's leading reciprocal file brands.

According to studies, T-Endo must files have been found to be more successful in torsional resistance and cyclic fatigue tests than many Reciprocal canal files.

**must™**  
multi use single technique

21 mm

25 mm

31 mm

tg*	13.04	
m25	25.06	
m40	40.04	
m50	50.04	

DC203TG	DC303TG	DC403TG
DC20325	DC30325	DC40325
DC20340	DC30340	DC40340
DC20350	DC30350	DC40350

must Asorti 4-file pack

Asorti	(1x TG, 1x M25, 1x M40, 1x M50)	DC203AST	DC303AST	DC403AST
--------	---------------------------------	----------	----------	----------

## Related GP



m25	DCGMUST-25
m40	DCGMUST-40
m50	DCGMUST-50

## Related GP



m25	DCPMUST-25
m40	DCPMUST-40
m50	DCPMUST-50



### Merkez

Molla Gürani Mah. Oğuzhan  
Cd.No: 45 Fatih/İstanbul / TÜRKİYE



### Fabrika

Akçaburgaz Mahallesi 3028.  
Sokak 9 Yapı Sanayi Sitesi No:1/33  
Esenyurt / İstanbul / Türkiye



[www.dentac.com](http://www.dentac.com)



[info@dentac.com](mailto:info@dentac.com)



+90 212 621 51 40

**dentac**<sup>®</sup>  
feel the innovation

