

Endless Love 

To Guard The New Life Safe



BM-100A

Transcutaneous Jaundice Detector



Stock code [300314]



Brief Introduction

■ Effective, User-friendly, Intelligent transmission

3.0 inch color touch screen, large character indicates and read measurement value clearly and rapidly;

On-screen display of multiple measurement data, and simultaneous storage of multiple patient information, achieving data management to patients.

■ Long endurance, environment friendly and durable

Equipped with charging base, enduring built-in chargeable battery, measuring 2000 times after battery is full;

Xenon flash lamp with long lifespan (measurement $\geq 150,000$ times);

Charging base with the function of integrated light checker, consumable-free, saving cost for check.


■ Accurate measurement and quick check

Effective check to identify patient with high-risk jaundice precisely, conform to the newest clinical guidance;

Constant measurement to multiple patients without repeating calibration;

Average measurement: this function permits 1~5 times setting of average measurement to ensure the accuracy on results;

Converting between two kinds of measurement unit mg/dL, $\mu\text{mol/L}$.



Standard Configuration

Consist of main body, base and power adapter.



Ningbo David Medical Device Co.,Ltd