



BIO STICK MICRO APPLICATORS

DISPOSABLE MICRO APPLICATORS



CONFIGURATION OVERVIEW

Material:	Flexible handle made of starch, tip is made of Nylon – please discard separately
Type:	Non-Sterile
Description:	<ul style="list-style-type: none"> - Ideal for application of sealants, bonding agents, etchants, hemostatic solutions, cavity liners and disclosing solutions. - Non-absorbent, non-linting fibers suspend solutions as small as 1/8 of a drop without spilling, dripping and waste. - Bendable neck "stays put" to facilitate the placement of material into difficult areas. - Flexible plastic handle for pinpoint accuracy application of any liquids

THECNICAL SPECIFICATION

Size	Regular blue	Fine yellow	Superfine beige
Length of the applicator (including the fiber tip) (mm)	100	100	100
Head Size (mm)	3.0	1.9	1.5
Cylinder (mm)	1.4	3.5	1.4

Conform to EUR REG MDR 745/2017 : Class 1 Medical Device according to rule 5 Annex VIII

Conform to ISO 10993-1/5/10, EN ISO 14971:2012, EN 15223-1:2016, EN 1041:2008+A1:2013

STORAGE AND HANDLING

Store in a well-ventilated room, temperature control $\leq 40^{\circ}\text{C}$, moisture $\leq 80\%$. Keep away from acid and avoid direct sunlight during storage and transportation. This product can be stored for two years without being disinfected and properly stored in a suitable environment.